

# Numerical Characterization of Ultrasound Elastography for the Early Detection of Deep Tissue Injuries

A thesis submitted in partial fulfilment of the requirements for the degree of Master of Science

by

Kenton David Hamaluik

# Contents

## 1. Introduction

The reasons for and goals of this research

## 2. Quasi-Static Ultrasound Elastography (QS USE)

Estimating stiffness using manual palpation

## 3. Acoustic Radiation Force Impulse (ARFI) Imaging

Using transducer-generated forces instead of manual palpation

## 4. Shear Wave Speed Quantification

Quantifying tissue stiffness using shear wave speeds


## 5. Conclusions

Recommendations and final thoughts

# Introduction

# Pressure Ulcers


- Pressure ulcers are secondary injuries
  - People with reduced mobility
- Skin breakdown due to moisture, shear / friction
- Categorized by NPUAP in stages
  - From shallow to deep



© National Pressure Ulcer Advisory Panel, used with permission.

# Deep Tissue Injuries

- Not all PU form “top-to-bottom”
  - Deep tissue injuries (DTI) form “bottom-to-top”
  - Eventually break out into stage III – IV pressure ulcers
- Tissue damage due to pressure and deformation
- Almost impossible to detect clinically



© National Pressure Ulcer Advisory Panel,  
used with permission.

# Deep Tissue Injury Detection

- ❖ T<sub>2</sub><sup>\*</sup>-weighted MRI in research settings
- ❖ Risk assessment scales in clinical settings
  - ❖ Norton, Braden, and Risk Assessment Pressure Sore scales

# Filling the Gaps



	DTI	B-Mode	QS USE	ARFI	Shear	FEM	Phantom	Animals	Humans	Characterization	Clinical
PU Risk scales	X	X	X	X	X	X	X	X	✓	X	✓
$T_2^*$ MRI	✓	—	—	—	—	✓	✓	✓	✓	X	X
Aoi et al.	✓	✓	X	X	X	X	X	X	✓	X	✓*
Deprez et al.	✓	X	✓	X	X	✓	✓	✓	X	X	✓
This work	✓	X	✓	✓	✓	✓	✓	X	X	✓	✓

# What?

# **Quasi-Static Ultrasound Elastography**


# Introduction

- Earliest form of ultrasound elastography
- Apply manual pressure to tissue
  - Measure localized deformation of tissue
- Magnitude of deformation related to stiffness
  - $\downarrow$  deformation  $\approx \uparrow$  stiffness  $\approx \uparrow$  damage magnitude




# Tracking Localized Deformation

- ✚ “Noise” isn’t actually noise
  - Scattering centres anchored in tissue
- ✚ Track motion of scattering centres between pre/post compression
  - Under assumptions of motion




Pre- and Post- Compression B-Mode Images of DTI

# Investigated Models



# Sample Resultant Elastogram



# Sample Characterization Results

- Show 4x4 of characterization plots

# Quasi-Static USE Outcomes

- ▣ QS USE **is** capable of DTI detection
- ▣ Detection sensitivity less than desirable
- ▣ Manual palpation is not ideal
  - ❖ Suggest **ARFI imaging** for machine-induced tissue deformation instead

# Acoustic Radiation Force Impulse Imaging

# Introduction

- QS USE has low detection sensitivity, manual palpation is difficult and unreliable
- ARFI imaging works on the same principles of QS USE
  - But uses transducer-generated force to displace tissue
- ↑ repeatability, ↑ inter-operator reliability
- By imparting acoustic energy to the tissue, body force is generated:

$$|\vec{F}| = \frac{2\alpha I}{c}$$

# How ARFI Imaging Works

- Normal ultrasound is only a couple periods long ( $\approx 2$  ms)
- ARFI imaging uses continuous beams ( $\approx 100$  ms)



- Acoustic energy is absorbed by tissue and causes deformation

# Sample ARFI Results

- 4x4 of sample ARFI results

# ARFI Imaging Outcomes

- ARFI more reliable than QS USE

# **Shear Wave Speed Quantification**


# Introduction

- QS USE and ARFI only provide qualitative measures of stiffness
- Measuring shear wave speed allows quantifiable calculation of stiffness
  - Can be used to accurately track over time and give absolute references
- Uses ARFI pulses to generate shear waves in tissue
  - Measure shear wave speed → calculate stiffness:

$$c_T = \sqrt{\frac{\mu}{\rho}} \rightarrow \mu = c_T^2 \rho$$

# Measuring Shear Wave Speed

1. Induce ARFI at desired depth
  - Shear waves radiate outward, like “ripples in a pond”
2. Monitor deformation along a line extending from the focal point using B-mode
  - Calculate speed of shear wave along this line



# Conclusions

# Experimental Validations

- ❖ Experiments carried out using a Siemens ACUSON S2000™ ultrasound machine with a Siemens 9L4 transducer
- ❖ Used a CIRS QA Phantom 049 soft tissue model with stiffness-varying spherical lesions
- ❖ Compared parametrically-identical experimental results to simulation
- ❖ Got this (→) characterization curve
- ❖ It works!

# Comparing Methods

# Recommendations

# **Additional Slides**

# Additional Slides

- ▶ backup 1

# backup1

herp