

Numerical Characterization of Ultrasound Elastography for the Early Detection of Deep Tissue Injuries

by


Kenton David Hamaluik

Supervisors: Dr. W. Moussa and Dr. M. Ferguson-Pell
Examiner: Dr. R. Burrell
Chair: Dr. M. Zuo


Contents

Introduction


Contents



Contents




Contents



Introduction

Deep Tissue Injuries (DTI)

- ✚ Secondary injuries (limited mobility)
- ✚ Stage III – IV pressure ulcers
- ✚ Caused by pressure & deformation
- ✚ No clinical detection
- ✚ Severe health burdens
[Russo et al., 2008]
 - 500,000 – 2,000,000 annually (USA)
 - Increased mortality, length of stay




© National Pressure Ulcer Advisory Panel, used with permission.

Current DTI Detection

- Clinical settings: risk assessment scales
 - ▶ Norton, Braden, and Risk Assessment Pressure Sore scales
 - ▶ No actual detection, only risk assessment
- Research settings: MRI
 - ▶ Not clinically feasible
- New research:
 - ▶ B-mode imaging [Aoi et al., 2009]
 - ▶ Blood / urine markers [Makhsous et al., 2010]
 - ▶ Ultrasound elastography [Deprez et al., 2011]

Ultrasound Elastography

- Mechanical **stiffness changes** with injury formation and progression



Adapted from [Gefen, 2009], used with permission

- Ultrasound elastography is a technology which measures tissue **stiffness**

Literature Review

- The purpose of this research was to gain understanding of and characterize the use of ultrasound elastography toward DTI detection


	DTI	B-Mode US	QS USE	ARFI	Shear	FE Models	Gel Phantoms	Animals	Humans	Characterized	Clinical
PU Risk scales	X	—	—	—	—	X	X	X	✓	✓	✓
MRI	✓	—	—	—	—	✓	✓	✓	✓	X	X
Aoi, 2009	✓	✓	X	X	X	X	X	X	✓	X	?
Makhsous, 2010	✓	—	—	—	—	—	—	✓	X	X	✓
Deprez, 2011	✓	X	✓	X	X	✓	✓	✓	X	X	✓
→ This work	✓	X	✓	✓	✓	✓	✓	?	?	✓	✓

Quasi-Static Ultrasound Elastography*

* Accepted for publication as: K. Hamaluik, W. Moussa,
and M. Ferguson-Pell, “Numerical Characterization of
Quasi-Static Ultrasound Elastography for the Detection of
Deep Tissue Injuries.” *IEEE Transactions on Medical Imaging*




Theory of Operation

- Manual deformation
- Magnitude of deformation related to stiffness
 - ↓ deformation
≈ ↑ stiffness
≈ ↑ damage magnitude




Theory of Operation

- ✚ Manual deformation
- ✚ Magnitude of deformation related to stiffness
 - ↓ deformation
 - ≈ ↑ stiffness
 - ≈ ↑ damage magnitude
- ✚ Track scattering centres
 - Assumptions of motion [Brusseau et al., 2008]



Investigated Models



d : lesion depth

h : lesion altitude

ϕS : lesion size

f : interrogation frequency

ε_{app} : applied strain


δ_{sep} : separation distance

b_r : blur radius


b_ρ : cluster density

r_{bl} : cluster radii


Sample Resultant Elastogram




Sample Quasi-Static Results




Sample Quasi-Static Results



Sample Quasi-Static Results



Sample Quasi-Static Results



Quasi-Static USE Outcomes

- ✚ Small lesions undetectable
 - ✖ Larger → more accurate
- ✚ Not reliant on:
 - ✖ Depth
 - ✖ Interrogation frequency
- ✚ Apply $\lesssim 5\%$ strain
- ✚ Complicated geometry can affect results
- ✚ Manual palpation not reliable
 - ✖ ARFI imaging as an alternative
- ✚ Can detect DTI
- ✚ Not ideally sensitive

Acoustic Radiation Force Impulse Imaging

Introduction



- QS USE has low detection sensitivity, unreliable
- ARFI imaging similar to QS USE
 - Transducer-generated force instead of manual
- ↑ repeatability, ↑ inter-operator reliability
- Absorbed energy generates force:

$$|\vec{F}| = \frac{2\alpha l}{c}$$

\vec{F} : ARFI force
 α : Absorption coefficient
 l : Intensity
 c : Longitudinal wave speed

How ARFI Imaging Works


- Normal ultrasound is only a couple periods long (≈ 2 ms)
- ARFI imaging uses continuous beams (≈ 100 ms)




- Acoustic energy is absorbed by tissue and causes deformation

Simulating ARFI Loads



- Use k-space pseudo-spectral method to solve the wave equation
 - k-Wave MATLAB[®] toolbox




Sample ARFI Results




Sample ARFI Results



Sample ARFI Results



Sample ARFI Results



ARFI Imaging Outcomes

- ARFI sensitivity \approx QS USE sensitivity
- ARFI not as dependent on lesion size as QS USE
- Not reliant on depth
- Complicated geometry can affect results
- Significant safety considerations
- ARFI more reliable than QS USE
- ARFI has limited penetration depth

Shear Wave Speed Quantification

Introduction


- QS USE and ARFI provide **qualitative** measures
- Shear wave speed quantification is **quantitative**
 - Absolute values
 - Track over time
- Uses ARFI pulses to generate shear waves
 - Measure shear wave speed → calculate stiffness:

$$\mu = c_T^2 \rho$$

μ : Shear modulus
 c_T : Shear wave speed
 ρ : Tissue density


Measuring Shear Wave Speed

1. Induce ARFI at desired depth
 - Shear waves radiate outward, like “ripples in a pond”




Measuring Shear Wave Speed

1. Induce ARFI at desired depth
 - Shear waves radiate outward, like “ripples in a pond”
2. Monitor deformation along a line extending from the focal point using B-mode ultrasound




Measuring Shear Wave Speed


1. Induce ARFI at desired depth
 - ❖ Shear waves radiate outward, like “ripples in a pond”
2. Monitor deformation along a line extending from the focal point using B-mode ultrasound
 - ❖ Calculate speed of shear wave along this line




Sample Shear Results




Sample Shear Results



Sample Shear Results



Sample Shear Results




Shear Wave Speed Outcomes

- ✚ Similar reliability to ARFI
- ✚ Difficult to detect:
 - Small lesions
 - Deep lesions
- ✚ Localized measures only
 - QS USE and ARFI provide domain interrogation
- ✚ More sensitive than ARFI and QS USE
- ✚ Quantifiable
 - Can monitor injury over time


Experimental Validation & Method Comparisons

Experimental Setup





Experimental Validations

- Siemens ACUSON S2000™ ultrasound machine with a Siemens 9L4 transducer and a CIRS QA Phantom 049 gel phantom
- Compared parametrically-identical experimental results to simulation



Comparing Methods



Conclusions

Recommendations

- ✚ Shear wave speed quantification yields best results
 - ❖ Use when possible
- ✚ ARFI more reliable than QS USE
 - ❖ Lacks penetration depth
- ✚ Small lesions ($r \lesssim 0.5$ cm) difficult to detect
 - ❖ MRI'd lesions: $r \gtrsim 1.0$ cm
- ✚ Future work: animal and human studies
 - ❖ Toward eventual clinical implementation

Conclusions

- DTI are a major health concern
 - ▶ No early detection
- Ultrasound elastography **can** detect DTI
 - ▶ Within limits
 - ▶ Clinically feasible
- Numerically characterized 3 modalities
 - ▶ Quasi-Static Ultrasound Elastography
 - ▶ Acoustic Radiation Force Impulse Imaging
 - ▶ Shear Wave Speed Quantification ←

Thank You

- ❖ Supervisors:
 - ❖ Dr. W. Moussa
 - ❖ Dr. M. Ferguson-Pell
- ❖ Lab-mates



Additional Slides

- ▣ References
- ▣ Derivation of ARFI
- ▣ Quasi-Static Parameters
- ▣ ARFI Parameters
- ▣ Shear Parameters
- ▣ Phantom Properties
- ▣ Quasi-Static Finite-Element Details
- ▣ Quasi-Static Mesh Dependency
- ▣ ARFI / Shear Finite-Element Details
- ▣ ARFI / Shear Mesh Dependency

References

- [Aoi et al., 2009] Aoi, N., Yoshimura, K., Kadono, T., Nakagami, G., Iizuka, S., Higashino, T., Araki, J., Koshima, I., and Sanada, H. (2009). **Ultrasound assessment of deep tissue injury in pressure ulcers: possible prediction of pressure ulcer progression.** *Plastic and reconstructive surgery*, 124(2):540–550.
- [Brusseau et al., 2008] Brusseau, E., Kybic, J., Deprez, J.-F. F., and Basset, O. (2008). **2-D locally regularized tissue strain estimation from radio-frequency ultrasound images: theoretical developments and results on experimental data.** *IEEE transactions on medical imaging*, 27(2):145–160.
- [Deprez et al., 2011] Deprez, J.-F. F., Brusseau, E., Fromageau, J., Cloutier, G., and Basset, O. (2011). **On the potential of ultrasound elastography for pressure ulcer early detection.** *Medical physics*, 38(4):1943–1950.
- [Gefen, 2009] Gefen, A. (2009). **Deep tissue injury from a bioengineering point of view.** *Ostomy/wound management*, 55(4):26–36.
- [Makhsous et al., 2010] Makhsous, M., Lin, F., Pandya, A., Pandya, M. S., and Chadwick, C. C. (2010). **Elevation in the serum and urine concentration of injury-related molecules after the formation of deep tissue injury in a rat spinal cord injury pressure ulcer model.** *PM & R : the journal of injury, function, and rehabilitation*, 2(11):1063–1065.
- [Russo et al., 2008] Russo, C., Steiner, C., and Spector, W. (2008). **Hospitalizations related to pressure ulcers among adults 18 years and older, 2006: Statistical brief no. 64.** *Agency for Health Care Policy and Research (US)*.

Derivation of Acoustic Force

$$\sigma_{ij,j} + \rho b_i = \rho f_i$$

$$\vec{F} = \nabla p_2 - \mu_{tissue} \nabla^2 \vec{v}_2$$

$$\vec{F} = \rho (\vec{v}_1 \nabla \cdot \vec{v}_1 + \vec{v}_1 \nabla \cdot \vec{v}_1)$$

$$\vec{F} = 2\rho \langle \vec{v} \vec{v},_x \rangle$$

$$\vec{v} = i\omega A e^{-\alpha x + i(\omega t - kx)} \hat{x}$$

$$|\vec{F}| = A^2 e^{-2\alpha x} \rho \alpha$$

$$\boxed{|\vec{F}| = \frac{2\alpha I}{c}}$$

- $\langle \rangle$: Time-average
- \vec{F} : Acoustic Force
- \vec{v} : Particle velocity
- ρ : Tissue density
- α : Absorption coefficient
- I : Acoustic intensity
- c : Speed of sound

QS USE Parameters

QS USE Parametric Study Parameters

Parameter	Symbol	Values	Units
Lesion depth	d	3.5, 6.5, 8.5 and 10.0	cm
Lesion altitude	h	1.25, 2.50 and 3.75	cm
Lesion diameter	ϕS	0.5, 1.0, 2.0 and 2.5	cm
Lesion stiffness ratio	E_{rel}	0.32, 0.56, 1.80 and 3.20	–
Ultrasound frequency	f	2, 4 and 8	MHz
Transducer-applied strain	ε_{app}	2.5, 5.0 and 10.0	%
Co-located separation distance	δ_{sep}	1.25, 1.50, 1.75 and 2.00	cm
Blurred lesion blur radius	b_r	1.0, 2.5, 5.0 and 7.5	mm
Clustered lesion density	b_p	10, 20, 30 and 40	cm^{-2}
Clustered lesion radius	r_{bl}	0.5, 1.0 and 1.5	mm
Visible human lesion width	ϕL	0.5, 1.0, 2.0 and 2.5	cm
Visible human lesion depth	d	6.25, 6.75 and 7.25	cm

ARFI Parameters

K-Space Psuedospectral Model Parameters

Property	Symbol	Value	Units
Nonlinearity parameter	B/A	8	-
Power law prefactor	α_0	0.7	$Np \text{ (rad/s)}^{-y} \text{ m}^{-1}$
Power law exponent	y	0.95	-
Density	ρ_0	1060	kg m^{-3}

ARFI Parameters

Simulated Material Parameters

Property	Symbol	Value	Units
Bulk Modulus	K	515.7	kPa
Shear Modulus	μ_{tissue}	1.0	kPa
Density	ρ	1060	kg m^{-3}

Generalized Maxwell Viscoelastic Material Model Parameters

Branch	Shear Modulus (Pa)	Relaxation Time (s)
1	791.0	2
2	66.5	40
3	0.6	80

ARFI Parameters

ARFI Parametric Study Parameters


Parameter	Symbol	Values	Units
ARFI interrogation frequency	f	1, 2, 4 and 6	MHz
Transducer width	w_{trans}	4, 8 and 10	cm
ARFI pulse cycles	n_c	3, 100, 300, 500 and 700	-
ARFI source pressure	P_{source}	4, 5, 6, 7 and 8	MPa
Lesion depth	d	1, 2, 3, 4, 5, 6, 7, 8 and 9	cm
Lesion diameter	ϕS	0.5, 1.0, 2.0 and 2.5	cm
Lesion stiffness ratio	E_{rel}	0.32, 0.56, 1.80 and 3.20	-
Blurred lesion blur radius	b_r	1.0, 2.5, 5.0 and 7.5	mm
Clustered lesion density	b_p	10, 20, 30 and 40	cm^{-2}
Clustered lesion radius	r_{bl}	0.5, 1.0 and 1.5	mm
Visible human lesion width	ϕL	0.5, 1.0, 2.0 and 2.5	cm

Shear Parameters

Shear Parametric Study Parameters

Parameter	Symbol	Values	Units
Lesion depth	d	1, 2, 3, 4, 5, 6, 7, 8 and 9	cm
Lesion diameter	ϕS	0.5, 1.0, 2.0 and 2.5	cm
Lesion offset	d_{off}	0.00, 1.25, 2.50 and 3.75	cm
Lesion stiffness ratio	E_{rel}	0.32, 0.56, 1.80 and 3.20	
Blurred lesion blur radius	b_r	1.0, 2.5, 5.0 and 7.5	mm
Clustered lesion density	b_ρ	10, 20, 30 and 40	cm^{-2}
Clustered lesion radius	r_{bl}	0.5, 1.0 and 1.5	mm
Visible human lesion width	ϕL	0.5, 1.0, 2.0 and 2.5	cm

CIRS Phantom Properties




CIRS Phantom Material Properties and Geometry

Property	Symbol	Value	Units
Nominal basal elastic modulus	E_{tissue}	25	kPa
Lesion elastic modulii	E_{tissue}		
Type I		4	kPa
Type II		14	kPa
Type III		45	kPa
Type IV		80	kPa
Speed of sound	c_0	1540	m s^{-1}
Acoustic attenuation	α	0.5	$\text{dB cm}^{-1} \text{MHz}^{-1}$


Quasi-Static FEA Details

- ✚ Quadratic Lagrangian triangular elements
- ✚ 9,468 degrees of freedom
- ✚ Solved using COMSOL Multiphysics




The Quasi-Static Mesh

Quasi-Static Mesh Dependency



ARFI / Shear FEA Details

- Quadratic Lagrangian triangular elements
- 51,096 degrees of freedom
- 25 ns constant time stepping
- Simulation stopped after focal displacement reached 1 % of maximal value
- Solved using COMSOL Multiphysics



The ARFI / Shear Mesh

ARFI / Shear Mesh Dependency

