

# Numerical Characterization of Ultrasound Elastography for the Early Detection of Deep Tissue Injuries

A thesis submitted in partial fulfilment of the requirements for the degree of Master of Science

by

Kenton David Hamaluik

# Contents

## 1. Introduction

The reasons for and goals of this research

## 2. Quasi-Static Ultrasound Elastography (QS USE)

Estimating stiffness using manual palpation

## 3. Acoustic Radiation Force Impulse (ARFI) Imaging

Using transducer-generated forces instead of manual palpation

## 4. Shear Wave Speed Quantification

Quantifying tissue stiffness using shear wave speeds


## 5. Conclusions

Recommendations and final thoughts

# Introduction

# Deep Tissue Injuries

- Deep tissue injuries (DTI)
  - Secondary injuries for those with limited mobility
  - Form deep in tissue
  - Eventually break out into stage III – IV pressure ulcers
- Tissue damage due to **pressure** and **deformation**
- Almost impossible to detect clinically
- Severe health and monetary burdens




© National Pressure Ulcer Advisory Panel,  
used with permission.

# Current DTI Detection

- Clinical settings: risk assessment scales
  - Norton, Braden, and Risk Assessment Pressure Sore scales
  - No actual detection, only risk assessment
- Research settings:  $T_2^*$ -weighted MRI
  - Not feasible clinically
- New research:
  - B-mode imaging (Aoi et al. [1])
  - Blood / urine markers (Mahksous et al. [2])
  - Ultrasound elastography (Deprez et al. [3])

# Ultrasound Elastography

- Mechanical **stiffness changes** with injury formation and progression



Adapted from Gefen [4], used with permission

- Ultrasound elastography is a technology which measures tissue **stiffness**

# Literature Review


	DTI	B-Mode US	QS USE	ARFI	Shear	FE Models	Gel Phantoms	Animals	Humans	Characterized Clinical
PU Risk scales	X	—	—	—	—	X	X	X	✓	✓
$T_2^*$ MRI	✓	—	—	—	—	✓	✓	✓	✓	X X
Aoi et al. [1]	✓	✓	X	X	X	X	X	X	✓	X X
Deprez et al. [3]	✓	X	✓	X	X	✓	✓	✓	X	X ✓
Mahksous et al. [2]	✓	—	—	—	—	—	—	✓	X	X ✓
This work	✓	X	✓	✓	✓	✓	✓	X	X	✓ ✓

- The purpose of this research was to gain understanding of and characterize the use of ultrasound elastography toward DTI

# **Quasi-Static Ultrasound Elastography**


# Introduction

- Earliest form of ultrasound elastography
- Apply manual pressure to tissue
  - Measure localized deformation of tissue
- Magnitude of deformation related to stiffness
  - $\downarrow$  deformation  $\approx \uparrow$  stiffness  $\approx \uparrow$  damage magnitude




# Tracking Localized Deformation

- + “Noise” isn’t actually noise
  - Scattering centres anchored in tissue
- + Track motion of scattering centres between pre/post compression
  - Under assumptions of motion (Brusseau et al. [5])




Pre- and Post- Compression B-Mode Images of DTI


# Investigated Models



# Sample Resultant Elastogram



# Sample Quasi-Static Results



# Quasi-Static USE Outcomes

- QS USE **is** capable of DTI detection
- Detection sensitivity less than desirable
- Manual palpation is not ideal
  - Suggest **ARFI imaging** for machine-induced tissue deformation instead
- Small lesions almost undetectable
- Use  $\leq 5\%$  applied strain
- Depth, interrogation frequency don't affect detection ability

# Acoustic Radiation Force Impulse Imaging



# Introduction

- QS USE has low detection sensitivity, manual palpation is difficult and unreliable
- ARFI imaging works on the same principles as QS USE
  - But uses transducer-generated force to displace tissue
- ↑ repeatability, ↑ inter-operator reliability
- By imparting acoustic energy to the tissue, body force is generated:

$$|\vec{F}| = \frac{2\alpha I}{c}$$

# How ARFI Imaging Works

- Normal ultrasound is only a couple periods long ( $\approx 2$  ms)
- ARFI imaging uses continuous beams ( $\approx 100$  ms)



- Acoustic energy is absorbed by tissue and causes deformation


# Simulating ARFI Loads

- Use k-space pseudo-spectral method to solve the wave equation
  - k-Wave MATLAB® toolbox



$$\frac{\partial \vec{v}}{\partial t} = -\frac{1}{\rho_0} \nabla p$$

$$\frac{\partial p}{\partial t} = -(2\rho + \rho_0) \nabla \cdot \vec{v} - \vec{v} \cdot \nabla \rho_0$$

$$p = c_0^2 \left( \rho + \vec{u} \cdot \nabla \rho_0 + \frac{B}{2A} \frac{\rho^2}{\rho_0} - L\rho \right)$$



# Sample ARFI Results



# ARFI Imaging Outcomes

- ✚ ARFI more reliable than QS USE
  - ❖ Computer-controlled deformation force
- ✚ Safety implications
  - ❖ High intensity acoustic waves can damage tissue
  - ❖ Difficult to get deep penetration without damage
- ✚ No effect on detection sensitivity of lesion size for radii  $\geq$  2.5 mm
- ✚ Lesion depth only a factor in generating radiation force
  - ❖ Soft tissue absorption makes deep penetration difficult
- ✚ Better at interrogating complicated geometry than QS USE

# **Shear Wave Speed Quantification**


# Introduction

- ❖ QS USE and ARFI only provide **qualitative** measures of stiffness
- ❖ Measuring shear wave speed allows quantifiable calculation of stiffness
  - ❖ Can be used to accurately track over time and give absolute references
- ❖ Uses ARFI pulses to generate shear waves in tissue
  - ❖ Measure shear wave speed → calculate stiffness:

$$\mu = c_T^2 \rho$$


# Measuring Shear Wave Speed

1. Induce ARFI at desired depth
  - Shear waves radiate outward, like “ripples in a pond”




# Measuring Shear Wave Speed

1. Induce ARFI at desired depth
  - Shear waves radiate outward, like “ripples in a pond”
2. Monitor deformation along a line extending from the focal point using B-mode ultrasound




# Measuring Shear Wave Speed

1. Induce ARFI at desired depth
  - Shear waves radiate outward, like “ripples in a pond”
2. Monitor deformation along a line extending from the focal point using B-mode ultrasound




# Measuring Shear Wave Speed

1. Induce ARFI at desired depth
  - ▶ Shear waves radiate outward, like “ripples in a pond”
2. Monitor deformation along a line extending from the focal point using B-mode ultrasound
  - ▶ Calculate speed of shear wave along this line



# Sample Shear Results




# Shear Wave Speed Outcomes

- ✚ More accurate than ARFI and QS USE
- ✚ Quantifiable
  - ❖ Can monitor injury over time
- ✚ Localized interrogation instead of full domain imaging
- ✚ Necessary to locate ARFI focal point 1.25 cm – 2.50 cm away from the region of interest
- ✚ Lesions with radii  $\lesssim$  5 mm difficult to detect
- ✚ Depths  $\lesssim$  8 cm
- ✚ Not strongly dependent on complicated geometry



# Conclusions

# Experimental Validations

- Siemens ACUSON S2000™ ultrasound machine with a Siemens 9L4 transducer and a CIRS QA Phantom 049 gel phantom
- Compared parametrically-identical experimental results to simulation



# Comparing Methods



# Recommendations

- ❖ Shear wave speed quantification yields best results
  - ❖ Quantifiable
  - ❖ Most accurate
- ❖ ARFI and shear wave speed quantification are depth-limited (tissue safety)
  - ❖ Quasi-static elastography for deep ( $\gtrsim 5$  cm) regions
- ❖ Small lesions ( $r \lesssim 5$  mm) are difficult to detect
- ❖ Complicated geometries can affect results
- ❖ Future work should involve animal and human studies
  - ❖ This is the first stage in employing ultrasound elastography as a clinical deep tissue injury detection mechanism

# Thank-You

Thank-you to my supervisors, Dr. Waled Moussa and Dr. Martin Ferguson-Pell for their guidance, my labmates for support, and the Project SMART team for funding and support.



# Additional Slides

- ▶ References
- ▶ Experimental Setup
- ▶ Derivation of ARFI
- ▶ Quasi-Static Parameters
- ▶ ARFI Parameters
- ▶ Shear Parameters
- ▶ Phantom Properties


# References

- [1] N. Aoi, K. Yoshimura, T. Kadono, G. Nakagami, S. Iizuka, T. Higashino, J. Araki, I. Koshima, and H. Sanada, "Ultrasound assessment of deep tissue injury in pressure ulcers: possible prediction of pressure ulcer progression." *Plastic and reconstructive surgery*, vol. 124, no. 2, pp. 540–550, Aug. 2009. [Online]. Available: <http://dx.doi.org/10.1097/prs.0b013e3181addb33>
- [2] M. Makhsous, F. Lin, A. Pandya, M. S. Pandya, and C. C. Chadwick, "Elevation in the serum and urine concentration of injury-related molecules after the formation of deep tissue injury in a rat spinal cord injury pressure ulcer model." *PM & R : the journal of injury, function, and rehabilitation*, vol. 2, no. 11, pp. 1063–1065, Nov. 2010. [Online]. Available: <http://view.ncbi.nlm.nih.gov/pubmed/21093844>

# References (cont.)

- [3] J.-F. F. Deprez, E. Brusseau, J. Fromageau, G. Cloutier, and O. Basset, "On the potential of ultrasound elastography for pressure ulcer early detection." *Medical physics*, vol. 38, no. 4, pp. 1943–1950, Apr. 2011. [Online]. Available: <http://view.ncbi.nlm.nih.gov/pubmed/21626927>
- [4] A. Gefen, "Deep tissue injury from a bioengineering point of view." *Ostomy/wound management*, vol. 55, no. 4, pp. 26–36, Apr. 2009. [Online]. Available: <http://view.ncbi.nlm.nih.gov/pubmed/19387094>
- [5] E. Brusseau, J. Kybic, J.-F. F. Deprez, and O. Basset, "2-D locally regularized tissue strain estimation from radio-frequency ultrasound images: theoretical developments and results on experimental data." *IEEE transactions on medical imaging*, vol. 27, no. 2, pp. 145–160, Feb. 2008. [Online]. Available: <http://dx.doi.org/10.1109/tmi.2007.897408>

# Experimental Setup



# Derivation of Acoustic Force

$$\sigma_{ij,j} + \rho b_i = \rho f_i$$

$$\vec{F} = \nabla p_2 - \mu_{tissue} \nabla^2 \vec{v}_2$$

$$\vec{F} = \rho \langle \vec{v}_1 \nabla \cdot \vec{v}_1 + \vec{v}_1 \nabla \cdot \vec{v}_1 \rangle$$

$$\vec{F} = 2\rho \langle \vec{v} \vec{v}_x \rangle$$

$$\vec{v} = i\omega A e^{-\alpha x + i(\omega t - kx)} \hat{x}$$

$$|\vec{F}| = A^2 e^{-2\alpha x} \rho \alpha$$

$$|\vec{F}| = \frac{2\alpha I}{c}$$

- $\langle \rangle$ : Time-average
- $\vec{F}$ : Acoustic Force
- $\vec{v}$ : Particle velocity
- $\rho$ : Tissue density
- $\alpha$ : Absorption coefficient
- $I$ : Acoustic intensity
- $c$ : Speed of sound

# QS USE Parameters

Parameter	Symbol	Values	Units
Lesion depth	$d$	3.5, 6.5, 8.5 and 10.0	cm
Lesion altitude	$h$	1.25, 2.50 and 3.75	cm
Lesion diameter	$\phi S$	0.5, 1.0, 2.0 and 2.5	cm
Lesion stiffness ratio	$E_{rel}$	0.32, 0.56, 1.80 and 3.20	–
Ultrasound frequency	$f$	2, 4 and 8	MHz
Transducer-applied strain	$\varepsilon_{app}$	2.5, 5.0 and 10.0	%
Co-located separation distance	$\delta_{sep}$	1.25, 1.50, 1.75 and 2.00	cm
Blurred lesion blur radius	$b_r$	1.0, 2.5, 5.0 and 7.5	mm
Clustered lesion density	$b_\rho$	10, 20, 30 and 40	$\text{cm}^{-2}$
Clustered lesion radius	$r_{bl}$	0.5, 1.0 and 1.5	mm
Visible human lesion width	$\phi L$	0.5, 1.0, 2.0 and 2.5	cm
Visible human lesion depth	$d$	6.25, 6.75 and 7.25	cm

# ARFI Parameters

Property	Symbol	Value	Units
Nonlinearity parameter	$B/A$	8	-
Power law prefactor	$\alpha_0$	0.7	$Np \text{ (rad/s)}^{-y} \text{ m}^{-1}$
Power law exponent	$y$	0.95	-
Density	$\rho_0$	1060	$\text{kg m}^{-3}$

# ARFI Parameters

Property	Symbol	Value	Units
Bulk Modulus	$K$	515.7	kPa
Shear Modulus	$\mu_{tissue}$	1.0	kPa
Density	$\rho$	1060	$\text{kg m}^{-3}$

Branch	Shear Modulus (Pa)	Relaxation Time (s)
1	791.0	2
2	66.5	40
3	0.6	80

# ARFI Parameters

Parameter	Symbol	Values	Units
ARFI interrogation frequency	$f$	1, 2, 4 and 6	MHz
Transducer width	$w_{trans}$	4, 8 and 10	cm
ARFI pulse cycles	$n_c$	3, 100, 300, 500 and 700	-
ARFI source pressure	$P_{source}$	4, 5, 6, 7 and 8	MPa
Lesion depth	$d$	1, 2, 3, 4, 5, 6, 7, 8 and 9	cm
Lesion diameter	$\phi S$	0.5, 1.0, 2.0 and 2.5	cm
Lesion stiffness ratio	$E_{rel}$	0.32, 0.56, 1.80 and 3.20	-
Blurred lesion blur radius	$b_r$	1.0, 2.5, 5.0 and 7.5	mm
Clustered lesion density	$b_\rho$	10, 20, 30 and 40	$\text{cm}^{-2}$
Clustered lesion radius	$r_{bl}$	0.5, 1.0 and 1.5	mm
Visible human lesion width	$\phi L$	0.5, 1.0, 2.0 and 2.5	cm

# Shear Parameters

Parameter	Symbol	Values	Units
Lesion depth	$d$	1, 2, 3, 4, 5, 6, 7, 8 and 9	cm
Lesion diameter	$\phi S$	0.5, 1.0, 2.0 and 2.5	cm
Lesion offset	$d_{off}$	0.00, 1.25, 2.50 and 3.75	cm
Lesion stiffness ratio	$E_{rel}$	0.32, 0.56, 1.80 and 3.20	
Blurred lesion blur radius	$b_r$	1.0, 2.5, 5.0 and 7.5	mm
Clustered lesion density	$b_p$	10, 20, 30 and 40	$\text{cm}^{-2}$
Clustered lesion radius	$r_{bl}$	0.5, 1.0 and 1.5	mm
Visible human lesion width	$\phi L$	0.5, 1.0, 2.0 and 2.5	cm

# CIRS Phantom Properties

Property	Symbol	Value	Units
Nominal basal elastic modulus	$E_{tissue}$	25	kPa
Lesion elastic modulus	$E_{lesion}$	8, 14, 45 and 80	kPa
Speed of sound	$c_0$	1540	$\text{m s}^{-1}$
Acoustic attenuation	$\alpha$	0.5	$\text{dB cm}^{-1} \text{MHz}^{-1}$
Lesion diameter	$\phi S$	10 and 20	mm
Lesion depth	$d$	15 and 35	mm